




Novatech
new biotechnology for life

a bess group company

Products for Pulmonology



STERITALC® LARGE PARTICLE SIZE TALC FOR USE IN THE PLEURAL CAVITY



STERITALC® consists of talcum which is mined in France and is specifically processed for medical use (talcum pleurodesis).

STERITALC® is suited for all indications of pleurodesis. It is non-soluble and induces permanent pleurodesis. Compared with tetracyclines, talcum is more effective and less painful.

As a rule, for malignant indications 3 to 5 g are used, for treatment of spontaneous pneumothorax 2 g are sufficient in most cases.

A critical side effect of talcum pleurodesis can be ARDS (Acute Respiratory Distress Syndrome). A possible cause of ARDS may be the systemic dissemination of talcum. In some cases, after application in the pleural cavity, talcum was found in other organs (kidneys, spleen, liver), too. The literature assumes that there is a relation between the talcum particle size and the systemic dissemination of talcum: smaller talcum particles appear to disseminate more than larger ones¹. The clinical picture also shows the effect of different particle sizes: talcum with a mean particle size below 15 µm induced stronger systemic and pulmonary inflammation reactions than talcum with a mean particle size of 25 µm².

STERITALC®, produced by NOVATECH, is specifically calibrated to a mean particle size of 25 µm in order to avoid systemic dissemination. Animal³ and clinical studies² show the lesser systemic dissemination.

A multi-center study showed that STERITALC® with its calibrated particle size can be safely used for pleurodesis of malignant pleural effusions. None of more than 550 patients developed ARDS⁴. The authors recommend to use no other talcum.

Another cause of ARDS may be a sepsis due to unsterile talcum, or talcum containing endotoxines⁵. This, too, can be excluded when STERITALC® is used, because STERITALC® is free of endotoxines and comes sterile.

¹ Ferrer, CHEST 2002; 122: 1018-1027

² Maskell, Am .J. Respir. Crit. Care Med. 2004; 170: 377-382

³ Fraticelli, CHEST 2002; 122:1737-1741

⁴ Janssen, Lancet 2007; 369: 1535-1539

⁵ Antony, Eur. Respir. J. 2001; 18: 402-419

Please see page 6 for more professional references.

STERITALC® LARGE PARTICLE SIZE TALC FOR USE IN THE PLEURAL CAVITY



STERITALC® comes in three dosage forms allowing different ways of application:

STERITALC® Vial F2, F4

- For use as a slurry (to be mixed with physiological saline solution; Xylocain may be added)
- For poudrage application via a thoracoscope / trocar

STERITALC® PF3

- For direct poudrage application via a thoracoscope / trocar (nebulization by air)
- No limitations regarding storage and handling
- Self-contained system - comes in a set with 3 g STERITALC®, a 420 mm cannula as well as a balloon nebulizer.

STERITALC® Spray

- For direct poudrage application via a thoracoscope / trocar (nebulization by aerosol)

REF	STERITALC® Description	Quantity of medical talc	Items/box
16903	F2: Vial, 50 ml	2 g	4
16913	F4: Vial, 50 ml	4 g	4
16833	Spray, with cannula 440 mm	3 g	2
16863	PF3: Vial 10 ml with balloon and cannula 420 mm (Poudrage Kit)	3 g	2 kits
16983	Accessory for 16863: Vial, 10 ml	3 g	4
STERILE			

„Es gilt als gesichert, dass es beim französischen Luzenac-Talc (... Steritalc®, NOVATECH) zu keiner systemischen Talkumdissemination kommt.“

(“It is regarded secure that the French Luzenac talcum (... Steritalc®, NOVATECH) does not lead to a systemic talcum dissemination.”)

Schnyder / Tschopp: Behandlung des Pneumothorax mittels internistisch-thorakoskopischer Talkumpleurodesse. Der Pneumologe 2010. 7: 357-363

“The most important clinical implication of our study is that large-particle talc can safely be used for pleurodesis.

Other talc preparations should not be used for this indication.”

Janssen et al.: Safety of pleurodesis with talc poudrage in malignant pleural effusion: a prospective cohort study. Lancet 2007; 369: 1535-1539

EWS® SILICONE SPIGOTS FOR SEGMENTAL AND SUBSEGMENTAL BRONCHI

Endobronchial Watanabe Spigots (EWS®) were developed by NOVATECH in close cooperation with Dr. Y. Watanabe, Okayama, Japan. They are made of medical grade silicone (implantable for more than 29 days) dyed with barium-sulfate for excellent radiodiagnostic visibility. EWS® are tapered and have an anatomical design with studs on the outside avoiding migration. They come sterile, individually blister-packed. Three sizes are available.

Indications

- persistent (even after successful pleural drainage), inoperable pneumothorax,
- bronchopleural fistula (especially following thoracic surgery) with a continuous loss of air despite thoracic suction drain, in cases in which surgical intervention is not indicated,
- temporary treatment of haemoptysis of peripheral origin in expectation of a bronchial arterial embolisation or surgery.¹⁾

A study²⁾ performed in Japan with 63 patients including 40 cases of intractable pneumothorax, 12 cases of pyothorax with fistula and 7 cases of pulmonary fistula has shown that using EWS® is safer and has more permanent positive results than conventional methods. After determination of the affected bronchi with a balloon catheter (alternatively X-ray etc.), the EWS® are placed with a flexible bronchoscope and forceps guided by the working tube of the bronchoscope. EWS® were successfully placed in 96.7% of the cases. The loss of air was stopped or significantly reduced in 77.6% of the cases. No severe complications occurred.



Three sizes are available:

- S = Ø 5 mm
- M = Ø 6 mm
- L = Ø 7 mm

Subsequent therapy

- For certain patients bronchial occlusion with EWS® can be envisaged as the only treatment. When it fails or the result is imperfect, pleurodesis (for example using STERITALC®) or other surgery (if not contra-indicated) can be considered.
- The spigots can be removed after the patient's condition has improved and the thoracic drain is removed. If there are difficulties in removing EWS® for any reason, removal is not necessary.

REF	EWS® Endobronchial Watanabe Spigots, sterile	
01EWS12A	12 EWS®	3 x S, 6 x M, 3 x L
01EWS3S	6 EWS®	6 x S (Ø 5 mm)
01EWS3M	6 EWS®	6 x M (Ø 6 mm)
01EWS3L	6 EWS®	6 x L (Ø 7 mm)
STERILE		

1) H. Dutau et al. Endobronchial Embolization with a Silicone Spigot as a Temporary Treatment for Massive Hemoptysis. *Respiration* DOI:101159/000092954, published online April 21, 2006

2) Watanabe Y. et al. Bronchial Occlusion with Endobronchial Watanabe Spigot, *J Bronchol.*, 10, 4, 2003

SINGLE USE BOUTIN TROCARS FOR THE PLEURA

For pleural biopsies and punctions

NOVATECH has launched trocars for the pleura in cooperation with Professor C. Boutin (Marseilles, France) in 1998. He redesigned the existing trocar to improve the quality of sampling.

The single use Boutin trocars have a working length of 78 mm. They come sterile and are available in two diameters (2 and 3 mm).

The three-sided mandrel and the improved shape of the hook of the 3 mm trocar give this device a sampling quality higher than with the Abrams needle¹⁾.

- 2 mm: for pleural puncture only
- 3 mm: for pleural biopsy and puncture



Boutin Trocar 2 mm: for pleural puncture only, with a lateral orifice and two mandrels: sharp resp. blunt.



Boutin Trocar 3 mm: for pleural biopsy and puncture, with lateral hook and two mandrels: sharp resp. blunt.



Features

- sterile, ready for use
- convenient dispenser box
- single use, minimizing risk of contamination
- less traumatic due to diameters adapted to intended use
- equipped with a three-way stopcock

REF	Description	Items per box
58023	Disposable BOUTIN Trocar for pleural puncture 2 mm, with two mandrels	12
58033	Disposable BOUTIN Trocar for pleural biopsy and puncture, 3 mm, with two mandrels	12
STERILE		

¹⁾ Study by Prof. C. BOUTIN - Hôpital La Conception, Marseilles
Presentation: CHEST/Toronto/1998

STERITALC® PROFESSIONAL REFERENCES

Antony VB; Loddenkemper R; Astoul P; Boutin C; Goldstraw P; Hott J; Rodriguez Panadero F; Sahn SA: Management of malignant pleural effusions. *European Respiratory Journal* 18 (2001): 402-419

Bloom AI; Wilson MW; Kerlan RK Jr; Gordon RL; LaBerge JM: Talc pleurodesis through small-bore percutaneous tubes. *Cardiovascular and Interventional Radiology* 22 (1999): 433-436

Bresticker MA; Oba J; LoCicero J; Greene R: Optimal pleurodesis: a comparison study. *The Annals of Thoracic Surgery* 55 (1993): 364-367

Cardillo G; Facciolo F; Carbone L; Regal M; Corzani F; Ricci A; Di Martino M; Martelli M: Long-term follow-up of video-assisted talc pleurodesis in malignant recurrent pleural effusions. *European Journal of Cardio-Thoracic Surgery* 21 (2002): 302-305

Diacon AH; Wyser C; Bolliger CT; Tamm M; Pless M; Perruchoud AP; Solér M: Prospective Randomized Comparison of Thoracoscopic Talc Poudrage under Local Anesthesia versus Bleomycin Instillation for Pleurodesis in Malignant Pleural Effusions. *Am. J. Respir. Crit. Care Med.* 162(4) (2000): 1445-1449

Erickson KV; Yost M; Bynoe R; Almond C; Nottingham J: Primary treatment of malignant pleural effusions: video-assisted thoracoscopic surgery poudrage versus tube thoracostomy. *The American Surgeon* 68 (2002): 955-959

Ferrer J; Montes JF; Villarino MA; Light RW; García-Valero J: Influence of particle size on extrapleural talc dissemination after talc slurry pleurodesis. *Chest* 122 (2002): 1018-1027

Ferrer J; Villarino MA; Tura JM; Traveria A; Light RW: Talc preparations used for pleurodesis vary markedly from one preparation to another. *Chest* 119 (2001): 1901-1905

Fratelli A; Robaglia-Schlupp A; Riera H; Monjanel-Mouterde S; Cau P; Astoul P: Distribution of calibrated talc after intrapleural administration: an experimental study in rats. *Chest* 122 (2002): 1737-1741

Gillissen A; Kellner S: Talkumpleurodesese mittels internistischer Thorakoskopie beim malignen oder chronischen Pleuraerguss. *Der Pneumologe* 5 (2010): 336-342

Glazer M; Berkman N; Lafair JS; Kramer MR: Successful talc slurry pleurodesis in patients with nonmalignant pleural effusion. *Chest* 117 (2000): 1404-1409

Janssen JP; Collier G; Astoul P; Tassi GF; Noppen M; Rodriguez-Panadero F; Loddenkemper R; Herth FJF; Gasparini S; Marquette CH; Becke B; Froudarakis ME; Driesen P; Bolliger CT; Tschopp JM: Safety of pleurodesis with talc poudrage in malignant pleural effusion: a prospective cohort study. *Lancet* 369 (2007): 1535 - 1539

Kennedy L; Harley RA; Sahn SA; Strange C: Talc Slurry Pleurodesis: Pleural Fluid and Histologic Analysis. *Chest* 107 (1995): 1707-1712

Kennedy L; Sahn SA: Talc Pleurodesis for Treatment of Pneumothorax and Pleural Effusion. *Chest* 106 (1994): 1215-1222

Keller SM: Current and future therapy for malignant pleural effusion. *Chest* 103 (1993): 635-675

Mager HJ; Maesen B; Verzijlbergen F; Schramel F: Distribution of talc suspension during treatment of malignant pleural effusion with talc pleurodesis. *Lung Cancer* 36 (2002): 77-81

Marom EM; Patz EF Jr; Erasmus JJ; McAdams HP; Goodman PC; Herndon JE: Malignant pleural effusions: treatment with small-bore-catheter thoracostomy and talc pleurodesis. *Radiology* 210 (1999): 277-281

Maskell NA; Kee, YCG; Gleeson FV; Hedley EL; Pengelly G; Davies RJO: Randomized Trials Describing Lung Inflammation after Pleurodesis with Talc of Varying Particle Size. *American Journal of Respiratory and Critical Care Medicine* 170 (2004): 377-382

Milanez de Campos JR; Filho LO; de Campos Werebe E; Sette H Jr; Fernandez A; Filomeno LT; Jatene FB: Thoracoscopy talc poudrage, a 15 Year experience. *Chest* 119 (2001): 801-806

Mourad IA; Abdel Rahman AR; Aziz SA; Saber NM; Fouad FA: Pleurodesis as a Palliative Treatment of Advanced Lung Cancer with Malignant Pleural Effusion. *Journal of the Egyptian Nat. Cancer Inst.* 16(3) (2004): 188-194

Nasreen N; Mohammed KA; Dowling PA; Ward MJ; Galffy G; Antony VB: Talc induces apoptosis in human malignant mesothelioma cells in vitro. *American Journal of Respiratory and Critical Care Medicine* 161 (2000): 595-600

Sanchez-Armengol A; Rodriguez-Panadero F: Survival and talc pleurodesis in metastatic pleural carcinoma, revisited. *Chest* 104 (1993): 1482-1485

Schnyder JM; Tschopp JM: Behandlung des Pneumothorax mittels internistisch-thorakoskopischer Talkumpleurodese. *Der Pneumologe* 5 (2010): 357-363

Schramel FM; Sutedja TG; Braber JC; van Mourik JC; Postmus PE: Cost effectiveness of video-assisted thoracoscopic surgery versus conservative treatment for first time or recurrent spontaneous Pneumothorax. *European Respiratory Journal* 9 (1996): 1821-1825

Tschopp JM; Bolliger CT; Boutin C: Treatment of spontaneous pneumothorax: why not simple talc pleurodesis by medical thoracoscopy. *Respiration* 67 (2000): 108-111

Tschopp JM; Boutin C; Astoul P; Janssen JP; Grandin S; Bolliger CT; Delaunois L; Driesen P; Tassi G; Perruchoud AP: Talcage by medical thoracoscopy for primary spontaneous pneumothorax is more cost-effective than drainage: a randomised study. *European Respiratory Journal* 20 (2002): 1003-1009

Viallat JR; Rey F; Astoul P; Boutin C: Thoracoscopic talc poudrage pleurodesis for malignant effusions. *Chest* 110 (1996): 1387-1393

Weissberg D; Ben-Zeev I: Talc Pleurodesis. *Journal of Thoracic and Cardiovascular Surgery* 106 (1993): 689-695

The products in this catalog are **CE**-marked.



NOVATECH SA, La Ciotat, France



Please note that only the instructions for use enclosed with the respective product apply. Details in this catalogue about the use of products serve as a guide only and reflect the information available at the time of print. If necessary, please request a current version!

Your local contact:

