

# Vyntus<sup>TM</sup> series with pneumotach sensor technology; International Organization for Standardization (ISO) 26782 dynamic waveform test results

Joachim Schwertberger, Florian Lotter and Dr. Herman Groepenhoff, Vyairé Medical

## Introduction

---

International guidelines<sup>1</sup> recommend the validation of clinical diagnostic spirometry systems against a known and accepted standard. Spirometers are now required to meet International Organization for Standardization (ISO) 26782 standards, which specifies requirements for spirometers intended for the assessment of pulmonary function in humans weighing more than 10 kg.<sup>2</sup> Consistent with other analyses in the literature<sup>3-5</sup>, we used the 13 dynamic waveforms required by ISO 26782<sup>2</sup> to drive a computer-controlled mechanical syringe for the validation of Vyntus<sup>TM</sup> SPIRO, Vyntus<sup>TM</sup> PNEUMO and Vyntus<sup>TM</sup> IOS\* Pulmonary Function Systems (Vyairé Medical, Mettawa, IL, USA) equipped with a pneumotach sensor using the software platform SentrySuite<sup>TM</sup> 3.20 (Vyairé Medical, Mettawa, IL, USA).

## Method

*Volume-time waveforms for evaluation of forced expired volume at 1 second (FEV1), forced expired volume at 6 seconds (FEV6), and forced vital capacity (FVC)*

A production spirometry device was connected through a MicroGard<sup>TM</sup> filter and ultrasonic sensor adaptor to a pump system for testing, oriented as it would be for human testing. The 13 ATS waveforms were generated by discharging the syringe into the spirometry device three times under ambient conditions. All readings were recorded.

## Accuracy test

The average of the three tests under ambient conditions is compared with the targeted standard value in the following way:

- Deviation = average – standard
- Percentage deviation =  $100 \cdot (\text{average} - \text{standard}) / \text{standard}$

The accuracy validation limits for volumes, which include the waveform generator inaccuracy according to ISO 26782<sup>2</sup>, were set to  $\pm 3.5\%$  of reading or  $\pm 0.050\text{L}$ , whichever is greater. However, current 2019 ERS/ATS spirometry guidelines narrowed these limits to  $\pm 2.5\%$  of reading.<sup>1</sup> These limits include the allowable inaccuracy for the pump system.

## Repeatability test


The FEV<sub>1</sub>, FEV<sub>6</sub>, and FVC data from the accuracy test are used to derive the range of the three recordings:

- Range = maximum – minimum


The repeatability validation limits for the volume measured at ambient conditions are  $\pm 2.5\%$  or  $\pm 0.050\text{L}$ , whichever is greater. A repeatability error occurs if the range or percentage range exceeds these limits.

## Results

Accuracy testing results (mean) are summarized in Figure 1 A-I showing values compared to those set by ERS/ATS.<sup>1</sup> Repeatability testing results are summarized in Figure 2 A-I.



**Figure 1** Accuracy of FEV<sub>1</sub>, FVC and FEV<sub>6</sub>



**Figure 2** Repeatability of  $FEV_1$ , FVC and  $FEV_6$

## Conclusion

Testing of the Vyntus™ SPIRO, Vyntus™ PNEUMO and Vyntus™ IOS with pneumotach sensor, in combination with the SentrySuite™ software platform, shows all ISO waveform and ERS/ATS requirements have been met with no failure (100% passing). From these results we can conclude that all three Vyntus systems with a pneumotach sensor in combination with the SentrySuite software platform passed all ISO standard waveform requirements.

## REFERENCES

1. Standardization of Spirometry 2019 Update An Official American Thoracic Society and European Respiratory Society Technical Statement.
2. International Organization for Standardization. ISO 26782 Anaesthetic and Respiratory Equipment – Spirometer Intended for the Measurement of Time Forced Expired Volumes in Humans. International Organization for Standardization Geneva, Switzerland; 2009.
3. Hankinson JL, Gardner RM. Standard waveforms for spirometer testing. Am Rev Respir Dis. Aug 1982;126(2):362-4. doi:10.1164/arrd.1982.126.2.362.
4. Miller MR, Hankinson J, Brusasco V, et al. Standardisation of spirometry. Eur Respir J. Aug 2005;26(2):319-38. doi:10.1183/09031936.05.00034805.
5. Nelson SB, Gardner RM, Crapo RO, Jensen RL. Performance evaluation of contemporary spirometers. Chest. Feb 1990;97(2):288-97. doi:10.1378/chest.97.2.288.

\*Vyntus Pneumo and IOS are not cleared for FDA 510K and are not available in the US market.

---

26125 North Riverwoods Blvd, Mettawa, IL 60045, USA

[vyaire.com](http://vyaire.com)

### For global distribution.

©2022 Vyaire Medical, Inc. or one of its affiliates. All rights reserved. Vyaire, the Vyaire Medical logo and all other trademarks are trademarks or registered trademarks of Vyaire Medical, Inc. or one of its affiliates. Vyaire's Medical devices are class I & IIa according to Medical Devices Directive 93/42/EEC or Medical Device Regulation EU 2017/745 as indicated on each declaration of conformity. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labeling. VYR-GBL-2100297

