

Pediatric – 3200 Series

The 3200 Series Pediatric Safe-T-Tubes™ feature an angled extraluminal limb to facilitate suctioning and airway management. Shortened limbs can be ordered directly from our facility, supplied non-sterile. Please allow up to 4 working days for custom orders.

Available in sizes 6 mm, 7 mm, 8 mm and 9 mm (outside diameter) in either clear or radiopaque implant-grade silicone. Each Safe-T-Tube™ set includes the tube, two Plug/Ring sets, Instructions for Use, and a PATIENT-PAK™ (including Cleaning & Care Instructions and a storage container for Plug/Ring sets). Additional Plug/Ring sets are available. Refer to the chart below for corresponding product codes.

Intended Use

To maintain an adequate tracheal airway, as well as to provide support in the stenotic trachea that has been reconstituted or reconstructed

Indications:

- Acute laryngotracheal injuries
- Intrathoracic tracheal stenosis
- Reconstructed trachea
- Reconstituted trachea
- Tracheal stenosis when the cervical or thoracic trachea cannot be repaired
- Segmental resection and anastomosis
- Cervical trachea that cannot be repaired

Supplied

Sterile

For single use

1 Tube, 2 Plug/Ring Sets, 1 IFU, 1 PATIENT-PAK™ per carton



includes 2 Plug/Ring sets



Features

- Angled extraluminal limb
- Implant-grade silicone
- Non-adherent smooth surface
- Safety ring and groove system
- Available in clear or radiopaque material

Size	REF		Dimensions (mm)						REF
	Clear	Radiopaque	A	B	C	D	E	F	Plug/Ring Set
6	320006	32006R	6	94	46	42	45	6	321066
7	320007	32007R	7	95	46	43	45	6	321066
8	320008	32008R	8	96	46	42	45	8	321086
9	320009	32009R	9	97	46	43	45	8	321086

