

Thoracic – 5200 Series

The 5200 Series Thoracic Safe-T-Tube™ features an extra-long lower limb for thoracic applications. Shortened limbs can be ordered directly from our facility, supplied non-sterile. Please allow up to 4 working days for custom orders.

Available in sizes 10 mm, 11 mm, 12 mm, 13 mm, 14 mm, 15 mm and 16 mm (outside diameter) in either clear or radiopaque implant-grade silicone. Each Safe-T-Tube™ set is supplied with the tube, two Plug/Ring sets, Instructions for Use, and a PATIENT-PAK™ (including Cleaning & Care Instructions and a storage container for Plug/Ring sets). Also available are 15 mm adaptors and additional Plug/Ring sets. Refer to chart below for corresponding product codes.

Intended Use

To maintain an adequate tracheal airway, as well as to provide support in the stenotic trachea that has been reconstituted or reconstructed

Indications

- Acute laryngotracheal injuries
- Intrathoracic tracheal stenosis
- Reconstructed trachea
- Reconstituted trachea
- Tracheal stenosis when the cervical or thoracic trachea cannot be repaired
- Segmental resection and anastomosis
- Cervical trachea that cannot be repaired

Supplied

Sterile

For single use

1 Tube, 2 Plug/Ring Sets, 1 IFU, 1 PATIENT-PAK™ per carton



includes 2 Plug/Ring sets



Features

- Extra-long distal limb
- Implant-grade silicone
- Non-adherent smooth surface
- Safety ring and groove system
- Available in clear or radiopaque material

Size	REF		Dimensions (mm)						REF	
	Clear	Radiopaque	A	B	C	D	E	F	Plug/Ring Set	15 mm Adaptor
10	520010	52010R	10	113	31	74	38	8	321106	321107
11	520011	52011R	11	115	30	74	50	11	321126	321127
12	520012	52012R	12	112	29	72	50	11	321126	321127
13	520013	52013R	13	116	31	74	55	11	321126	321127
14	520014	52014R	14	118	32	75	58	11	321126	321127
15	520015	52015R	15	119	33	75	56	11	321126	321127
16	520016	52016R	16	119	33	75	65	11	321126	321127

