

# Montgomery® Safe-T-Tubes™

## Tapered – 7200 Series

The 7200 Series Tapered Safe-T-Tube™ is designed to assist with the management of tracheal, subglottic, and glottic stenoses. The superior end conforms to the contour of the glottis and subglottis, while the larger inferior portion simultaneously acts as a stent for the trachea. The tapered end is designed to project up beyond the true vocal cords to the level of the laryngeal ventricle. Shortened limbs can be ordered directly from our facility, supplied non-sterile. Please allow up to 4 working days for custom orders.

Available in sizes 8 mm/10 mm and 10 mm/13 mm (outside diameters) in either clear or radiopaque implant-grade silicone. Each Safe-T-Tube™ set includes the tube, two Plug/Ring sets, Instructions for Use, and a PATIENT-PAK™ (including Cleaning & Care Instructions and a storage container for Plug/Ring sets). Also available are 15 mm adaptors and additional Plug/Ring sets. Refer to the chart below for corresponding product codes.

### Intended Use

To maintain an adequate tracheal airway, as well as to provide support in the stenotic trachea that has been reconstituted or reconstructed

### Indications

- Acute laryngotracheal injuries
- Intrathoracic tracheal stenosis
- Reconstructed trachea
- Reconstituted trachea
- Tracheal stenosis when the cervical or thoracic trachea cannot be repaired
- Segmental resection and anastomosis
- Cervical trachea that cannot be repaired

### Supplied

Sterile

For single use

1 Tube, 2 Plug/Ring Sets, 1 IFU, 1 PATIENT-PAK™ per carton



includes 2 Plug/Ring sets



### Features

- Tapered upper limbs
- Implant-grade silicone
- Non-adherent smooth surface
- Safety ring and groove system
- Available in clear or radiopaque material

Size	REF		Dimensions (mm)							REF	
	Clear	Radiopaque	A	B	C	D	E	F	G	Plug/Ring Set	15 mm Adaptor
8/10	721008	71008R	8	163	77	78	69	8	10	321106	321107
10/13	721310	71310R	10	163	76	76	72	11	13	321126	321127

