


# INFLATION DEVICES



Analog Inflation Device

## THE GOLD STANDARD INFLATION DEVICE FOR INTERVENTIONAL BALLOON INFLATION

With over 10 million devices used worldwide, busy interventional labs can rely on the basixCOMPAK™ for their procedures. The basixCOMPAK offers 30 ATM pressure and 20 mL volume and can be packaged with a wide variety of angioplasty kits and procedure packs.



For additional kit configurations, please visit [merit.com/compak](http://merit.com/compak).

Catalog Number	Product Description	Units Per Box
IN4130	basixCOMPAK	5
IN4330	basixCOMPAK with Honor	5
IN4352	basixCOMPAK with Honor + Extension Tubing	5
IN4530	basixCOMPAK with AccessPLUS	5
IN4952	basixCOMPAK with FLO50	5

**Indications for Use:** The basixCOMPAK Inflation Device is used to inflate and deflate an angioplasty balloon or other interventional device and to measure the pressure within the balloon.