




Understand. Innovate. Deliver.™

BIG60®

INFLATION DEVICE





Merit Medical has been an Inflation Device Market Leader for more than 25 Years



Merit Works Closely with Health Professionals to Produce the Highest Quality Products



10M

Merit's Inflation Devices have been Used in more than 10 Million Interventions Worldwide



Merit is a Top Producer of Inflation Devices

BIG60 INFLATION DEVICE

Merit Medical's BIG60 Inflation Device is designed to rapidly inflate and completely deflate balloon dilators while accurately monitoring and displaying crucial inflation pressures. Our ergonomically designed 60 mL inflation device is a one-piece, easy-to-use device that pairs Merit Medical's market-leading inflation device experience with our stent technology.

Data on file

INFLATION DEVICE FEATURES

CONVENIENT ALL-IN-ONE DESIGN

Requires no assembly or cleaning before use and is completely disposable

LOCK RELEASE HANDLE

Allows for easy and accurate one-handed control of inflation and deflation

ACCURACY

With the locking mechanism engaged, pressure can be fine-tuned

ANGLED ANALOG PRESSURE GAUGE

Glow in the dark and provides accurate pressure readings and high visibility in low-light environments

BIG60 Inflation Device

60 mL Inflation Device with Angled Gauge

CATALOG NUMBER

ENDO-AN6012 Packaged 5 per box



Merit Medical Europe, Middle East & Africa

Amerikalaan 42, 6199 AE Maastricht-Airport, Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West, Galway, Ireland
+35 39 170 37 33

Merit Medical France

136 Ave. Joseph Kessel, 78960 Voisins Le Bretonneux, France
+33 (0)1 30 57 24 24

Merit Medical Systems, Inc.

1600 West Merit Parkway, South Jordan, UT, 84095 USA
+1 801 253 1600

Asia, Pacific Rim, Latin America, Caribbean & Canada

+1 801 208 4438

Free phone for specific country

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (French)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Sweden
020 792 445

Sweden
020 792 445

Switzerland
00800 9000 5060

UK
0800 973 115