

# Embosphere® Microspheres

# EmboGold® Microspheres

# CATHETER COMPATIBILITY CHART


**Merit Maestro®**  
Microcatheter



**Merit Pursue™**  
Microcatheter



**SwiftNINJA®**  
Steerable Microcatheter



Compatible Microsphere Size (µm)		50-100	40-120	100-300	300-500	500-700	700-900	900-1200
Catheter	Inner Diameter							
Catheter	0.038" (0.97 mm)							
Catheter	0.035" (0.89 mm)							
EmboCath® Plus Microcatheter	0.028" (0.71 mm)							
Merit Maestro® Microcatheter 2.9F	0.027" (0.69 mm)							
Merit Maestro Microcatheter 2.8F	0.024" (0.62 mm)							
SwiftNINJA® Steerable Microcatheter 2.4F	0.021" (0.54 mm)							
Merit Maestro Microcatheter 2.4F	0.020" (0.52 mm)							
Merit Pursue™ Microcatheter 2.0F	0.020" (0.50 mm)							
Merit Maestro Microcatheter 2.1F	0.018" (0.46 mm)							
Merit Pursue Microcatheter 1.7F	0.016" (0.40 mm)							

Above catheter compatibility information based on internal testing.

**Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.**



Understand. Innovate. Deliver.™

[merit.com](http://merit.com)

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733