
## Long Introducers



Merit's Long Introducers have a braided steel shaft for excellent torqueability. The tapered dilator and soft, atraumatic tip is designed for enhanced tissue crossing. The radiopaque sheath and dilator enables visualization and aids in proper placement.



The Long Introducers work in unison with Merit's Long Needles.

FND-019-00

## FEATURES & BENEFITS

Braided Steel Core	For excellent torqueability
Lubricious Coating, Tapered Dilator and Soft, Atraumatic Tip	Enhances tissue crossing
Snap-Fit Dilator	Locks in place
Radiopaque Sheath and Dilator	Aids in radiopacity for visualization and proper placement

## Long Introducers

For a complete list of part numbers, visit MeritOEM.com/long-introducers