



The Merit MAK-NV™ utilizes a small co-axial introducer with dilator and stiffening cannula and guide wire for placement of larger diameter guide wires in non-vascular procedures.

- Excellent transition for a smooth introduction
- Radiopaque marker tip clearly shows tip of introducer
- Superior tip design decreases buckling and introducer peel back during insertion
- Echo-enhanced needle increases visualization for ultrasound placement
- 4F introducer ID allows 0.018" guide wire to be left in place for safety during larger wire placement



MERIT MAK-NV™ MINI ACCESS KIT NON-VASCULAR

Catalog Number	Introducer	Dilator	Guide Wire 0.018" (0.46 mm) x 60 cm	Guide Wire Heavy Duty/ Double Ended 0.038" (0.97 mm) 3 mm J	Needle
MAK-NV-001	6F x 20 cm	4F	SS/PT	150 cm	21G x 15 cm Trocar
MAK-NV-002	6F x 20 cm	4F	NT/PT	150 cm	21G x 15 cm Trocar
MAK-NV-003	6F x 20 cm	4F	SS/PT	–	21G x 15 cm Trocar
MAK-NV-004	6F x 20 cm	4F	NT/PT	–	21G x 15 cm Trocar
MAK-NV-005*	–	–	–	–	21G x 15 cm Trocar
MAK-NV-006*	6F x 20 cm	4F	–	–	21G x 15 cm Trocar
MAK-NV-007*	6F x 20 cm	4F	–	–	–
MAK-NV-008	6F x 20 cm	4F	SS/PT	80 cm	21G x 15 cm Trocar
MAK-NV-009	6F x 20 cm	4F	NT/PT	80 cm	21G x 15 cm Trocar
MAK-NV-010	6F x 20 cm	4F	SS/PT	150 cm	21G x 15 cm Chiba
MAK-NV-011	6F x 20 cm	4F	NT/PT	150 cm	21G x 15 cm Chiba
MAK-NV-012	6F x 20 cm	4F	SS/PT	–	21G x 15 cm Chiba
MAK-NV-013	6F x 20 cm	4F	NT/PT	–	21G x 15 cm Chiba
MAK-NV-014*	–	–	–	–	21G x 15 cm Chiba
MAK-NV-015*	6F x 20 cm	4F	–	–	21G x 15 cm Chiba
MAK-NV-016	6F x 20 cm	4F	SS/PT	80 cm	21G x 15 cm Chiba
MAK-NV-017	6F x 20 cm	4F	NT/PT	80 cm	21G x 15 cm Chiba
MAK-NV-018	6F x 20 cm	4F	SS/PT	150 cm	21G x 20 cm Trocar
MAK-NV-019	6F x 20 cm	4F	NT/PT	150 cm	21G x 20 cm Trocar
MAK-NV-020*	–	–	–	–	21G x 20 cm Trocar
MAK-NV-021*	–	–	–	–	21G x 20 cm Chiba
MAK005-60	–	–	SS/PT	–	–
MAK001N-60	–	–	NT/PT	–	–

SS = Stainless Steel | PT = Platinum Tip | NT = Nitinol | Packaged 10 per box. | *Packaged 5 per box.