

Our ReSolve® Biliary Locking Drainage Catheter is designed for fast, unimpeded drainage and offers up to 64% more drainage area¹ than other drainage catheters, with multiple large holes and different size configurations.

- Atraumatic internal pigtail designed to reduce tissue irritation
- Depth markers and radiopaque marker band for fast and accurate positioning
- Two drainage hole patterns to accommodate multiple lesion locations and varying patient anatomy

RESOLVE® BILIARY LOCKING DRAINAGE CATHETER

| Catalog Number | French Size | Working Length | Drainage Location | Pigtail to Distal Shaft Drainage Hole | Drainage Holes | Securement Device |
|----------------|--|----------------|-------------------|---------------------------------------|----------------|-------------------|
| RBC-8-038 | 8.5 | 40 cm | 14 cm | 5 cm | 17 | – |
| RBC-10-038 | 10 | 40 cm | 14 cm | 5 cm | 17 | – |
| RBC-12-038 | 12 | 40 cm | 14 cm | 5 cm | 18 | – |
| RBC-14-038 | 14 | 40 cm | 14 cm | 5 cm | 18 | – |
| RBDC-8-038 | 8.5 | 40 cm | 11 cm | 2 cm | 17 | – |
| RBDC-10-038 | 10 | 40 cm | 11 cm | 2 cm | 17 | – |
| RBDC-12-038 | 12 | 40 cm | 11 cm | 2 cm | 18 | – |
| RBDC-14-038 | 14 | 40 cm | 11 cm | 2 cm | 18 | – |
| RBC-8-SFX | 8.5 | 40 cm | 14 cm | 5 cm | 17 | StayFIX® (small) |
| RBC-10-SFX | 10 | 40 cm | 14 cm | 5 cm | 17 | StayFIX (small) |
| RBC-12-SFX | 12 | 40 cm | 14 cm | 5 cm | 18 | StayFIX (small) |
| RBC-14-SFX | 14 | 40 cm | 14 cm | 5 cm | 18 | StayFIX (large) |
| RBDC-8-SFX | 8.5 | 40 cm | 11 cm | 2 cm | 17 | StayFIX (small) |
| RBDC-10-SFX | 10 | 40 cm | 11 cm | 2 cm | 17 | StayFIX (small) |
| RBDC-12-SFX | 12 | 40 cm | 11 cm | 2 cm | 18 | StayFIX (small) |
| RBDC-14-SFX | 14 | 40 cm | 11 cm | 2 cm | 18 | StayFIX (large) |
| RBC-8-REV | 8.5 | 40 cm | 14 cm | 5 cm | 17 | Revolution™ |
| RBC-10-REV | 10 | 40 cm | 14 cm | 5 cm | 17 | Revolution |
| RBC-12-REV | 12 | 40 cm | 14 cm | 5 cm | 18 | Revolution |
| RBC-14-REV | 14 | 40 cm | 14 cm | 5 cm | 18 | Revolution |
| RBDC-8-REV | 8.5 | 40 cm | 11 cm | 2 cm | 17 | Revolution |
| RBDC-10-REV | 10 | 40 cm | 11 cm | 2 cm | 17 | Revolution |
| RBDC-12-REV | 12 | 40 cm | 11 cm | 2 cm | 18 | Revolution |
| RBDC-14-REV | 14 | 40 cm | 11 cm | 2 cm | 18 | Revolution |
| RC-RTOOL | Repositioning Tool Packaged 15 per box | | | | | |

Packaged 5 per box. | Max guide wire: 0.038" (0.97 mm)

