



Designed for drainage of abscess and miscellaneous drainage procedures, the ReSolve® Non-Locking Drainage Catheter comes in a J or straight tip and is packaged with a metal stiffening cannula and a trocar stylet.

- Large drainage holes to maximize drainage
- Tapered tip and hydrophilic coating designed for a smooth introduction to reduce patient discomfort
- Color-coded hubs that match French sizing charts for easy identification
- Flexible polyurethane material to prevent kinking and to ensure excellent catheter rebound¹

RESOLVE® NON-LOCKING DRAINAGE CATHETER

Catalog Number	French Size	Working Length	Tip Configuration	Drainage Holes
RNL-6-038S	6.5	20 cm	Straight	4
RNL-7-038S	7.5	25 cm	Straight	4
RNL-8-038S	8.5	25 cm	Straight	4
RNL-10-038S	10	25 cm	Straight	4
RNL-12-038S	12	25 cm	Straight	5
RNL-6-038J	6.5	20 cm	J Curve	4
RNL-7-038J	7.5	25 cm	J Curve	4
RNL-8-038J	8.5	25 cm	J Curve	4
RNL-10-038J	10	25 cm	J Curve	4
RNL-12-038J	12	25 cm	J Curve	5

Packaged 5 per box | Max guide wire: 0.038" (0.97 mm)

